

## ONE YEAR LIMITED WARRANTY

Obusforme by Conair will (at our option) repair or replace your appliance free of charge for 12 months from the date of purchase if the appliance is defective in workmanship or materials. To obtain service under this warranty, return the defective product to the service center nearest you, together with your sales slip and \$7.50 for postage and handling. In the absence of a purchase receipt, the warranty period will be 12 months from the date of manufacture. **ANY IMPLIED WARRANTIES, OBLIGATIONS OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 12 MONTH DURATION OF THIS WARRANTY.**

**IN NO EVENT SHALL OBUSFORME BY CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER.**

This warranty gives you specific legal rights, and you may also have other rights which vary from province to province.

*Obus Forme® and The Body Care Experts™ are registered trademarks owned by Integral Orthopedics, Inc. used under license.*

**Distributed by:**  
**CONAIR CONSUMER PRODUCTS INC.**

### **SERVICE CENTRE:**

Conair Consumer Products Inc.  
100 Conair Parkway,  
Woodbridge, ON  
L4H 0L2

OBHP08C

For more information on any Conair Products call:  
1-800-472-7606

Visit our Web site: [www.conaircanada.ca](http://www.conaircanada.ca)

Consumer Call Centre E-mail: [Consumer\\_Canada@Conair.com](mailto:Consumer_Canada@Conair.com)

©2007 Conair Consumer Products Inc.  
100 Conair Parkway, Woodbridge, ON L4H 0L2

**IB-7784**

07CN10253

# **IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY**

**DANGER**– TO REDUCE THE RISK OF BURNS, ELECTRIC SHOCK AND FIRE, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS.

1. READ ALL INSTRUCTIONS BEFORE USING.
2. DO NOT USE WHILE SLEEPING.  
DO NOT USE ON AN INFANT.
3. THIS PAD IS NOT TO BE USED ON OR BY AN INVALID, A PARA-

6. DO NOT USE IN OXYGEN ENRICHED ENVIRONMENT OR NEAR EQUIPMENT WHICH STORES OR EMITS OXYGEN.
7. **DO NOT USE PINS OR OTHER METALLIC MEANS TO FASTEN THIS PAD IN PLACE.**
8. **DO NOT SIT ON OR AGAINST OR CRUSH PAD - AVOID SHARP FOLDS.**
9. **THIS POWER UNIT IS INTENDED TO BE CORRECTLY ORIENTATED IN A VERTICAL OR FLOOR MOUNT POSITION.**
10. **NEVER PULL THIS PAD BY THE SUPPLY CORD AND DO NOT USE**

## **CAUTIONS**

- This massager should not be used in the following cases unless directed by a physician: Tuberculosis, benign and malignant tumors, hemorrhages, inflammations of the skin, phlebitis and thrombosis, open and fresh wounds, bruises, broken skin or varicose veins or undetermined calf pain or joint pain.
- Consult your physician, should you have any therapeutic questions regarding its use.
- In case of pregnancy, diabetes, or illness, consult your physician before using massager.
- Do not use on any unexplained pain before consulting your physician.
- Massage should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your physician.
- Do not fall asleep while using this massager.
- Keep away from children. This unit is not a toy. Children should not use it or play with it.

## **MAINTENANCE INSTRUCTIONS**

- 1) Your massager is designed for household use, and is virtually maintenance-free. No lubrication is needed. To clean the head and body of the massager, use a damp cloth only. Do not use any strong chemical cleaners such as Benzine.
- 2) If the cord or massager is damaged or does not operate properly, return it to an authorized service center. Any required servicing should be performed by a Conair authorized service centre.
- 3) Never allow the power supply cord to be excessively pulled or twisted. Inspect the cord frequently for damage. Stop use immediately if damage is visible, or unit stops, or operates intermittently. Never wrap the cord around the appliance.

DOES NOT FUNCTION SATISFAC-  
TORILY, RETURN TO AUTHORIZED  
SERVICE CENTER AT ADDRESS  
GIVEN.

15. DO NOT USE THIS PAD WITH  
LINIMENT, SALVE, OR OINTMENT  
PREPARATIONS THAT CONTAIN  
HEAT PRODUCING INGREDIENTS.  
SKIN BURNS COULD RESULT.
16. USE ONLY SUPPLIED TRANS-  
FORMER WITH THIS PAD.
17. SAVE THESE INSTRUCTIONS.

# GET TO KNOW YOUR MASSAGING HEATING PAD

**King Size Pad**  
11.5" x 24"

**Self-adhesive Belt**  
for easy application

**UL Listed  
Transformer  
Included**

**On/Off Heat and  
Off/Hi/Low Massage  
Option**

**Massaging Heating Pad**  
for custom comfort to ease  
muscle aches and pains



# OBUSFORME®

THE BODYCARE EXPERTS™



## Massaging Heating Pad

### Model OBHP08C Instruction Booklet

For your safety and continued enjoyment of this products, always read the instruction book carefully before using.

LYZED PERSON, A SLEEPING OR UNCONSCIOUS PERSON, A PERSON WITH DIABETES, OR A PERSON WITH POOR BLOOD CIRCULATION. DO NOT USE ON AREAS OF INSENSITIVE SKIN.

4. BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING. CHECK SKIN UNDER PAD FREQUENTLY TO AVOID BURNING AND BLISTERING.
5. PLACE PAD ON TOP OF AND NOT UNDER THE PART OF THE BODY NEEDING HEAT - BURNS CAN RESULT REGARDLESS OF CONTROL SETTING.



## **THE CORD AS A HANDLE.**

11. LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING MAY DAMAGE CORD AND INTERNAL PARTS.
12. DISCARD THE PAD IF COVER SHOWS ANY SIGNS OF DETERIORATION (SUCH AS BLISTERING OR CRACKING).
13. USE THE SUPPLY TRANSFORMER ONLY ON A 120 VOLT AC CIRCUIT. UNPLUG WHEN NOT IN USE.
14. DO NOT TAMPER WITH THIS PAD IN ANY WAY. THERE ARE NO USER SERVICEABLE PARTS.  
IF FOR ANY REASON THIS PAD

## **RECOMMENDED MASSAGING TIMES**

- Generally, massage should last 10 - 20 minutes on any muscle of the body. Treatment may be applied once or twice daily.
- As people's metabolisms varies, so does their response to massage and you will soon adjust the time to suit your particular needs and personal comfort. If you have any doubts regarding treatment, we recommend you consult your physician.

## **TO OPERATE**

To operate the massaging heating pad, first set both the massaging and heat switches to off and plug the transformer into a standard 120v electrical outlet. Plug the power jack into the hand controller. Now turn on either the massaging feature by switching it to HI or LO or the heat feature to "ON". You can use these 2 features together or alone.

For maximum benefit, use should not exceed 20 minutes at one time.

## **TIMER FEATURE**

The massage motors will automatically shut-off after approximately 20 minutes. To restart, first set the massage switch back to OFF and then set the switch to LO or HI again. Please note the heating pad operates independently of the timer.